



VISULAS 532s

Compact laser workstation
for state-of-the-art retinal therapy



We make it visible.




// PHOTOCOAGULATION
MADE BY CARL ZEISS

When darkness led to enlightenment,
when spontaneous retinal damage led to
an effective therapy.

*Sparked by the solar eclipse in 1945
and driven by visionary partners*

Carl Zeiss has been committed to precision, innovation and passion for more than half a century in retinal photocoagulation. Ever since the first photocoagulation with sunlight was tested from the rooftop of a Hamburg eye clinic in 1949, dedicated Zeiss employees have sought intense dialogue with visionary ophthalmologists of their time. Dialogue that led, in 1957, to the first commercially available photocoagulator in the world. Simultaneously this marked the birth of contactless surgery on the human eye.

Carl Zeiss has been a strong driving force behind the focused optical therapy of retinal diseases. Committed to this tradition, the VISULAS 532s ranks seamlessly into a chain of innovative solutions with one sole objective: precise, effective and gentle therapy to preserve the eyesight and quality of life of patients.



VISULAS 532s


Mastering everyday life in the clinic.
Reliable. Competent. Superior.

Exceptionally precise beam control and treatment accuracy.

Extremely compact and reliable in daily clinical practice.

Outstandingly ergonomic and intuitive to operate.

The VISULAS® 532s is an impressive, powerful, diode-pumped solid-state laser for controlled and gentle photocoagulation of the retina.



3D illustration of a single pulse coagulation



Fundus image of a single pulse coagulation
near the optic nerve head

Effective and reproducible treatment

The VISULAS 532s is setting performance benchmarks. The efficient solid-state laser reveals its sustainable strength even in high-demand practices. The VISULAS 532s has sufficient power reserves for any treatment strategy – whether it is with short laser pulses in the millisecond range for a gentle grid treatment, or long pulses for an effective retinopexy or even angioma sclerotherapy. Its built-in thermoelectric cooling system ensures maximum temporal stability of the laser power and thus meets the basic prerequisite for reproducible clinical results.

A laser workstation par excellence

In combination with the ZEISS Laser Slit Lamp LSL 532s, the VISULAS 532s becomes a perfectly integrated laser workstation.

The electronic micromanipulator allows for sensitive, synchronous laser beam positioning and slit lamp illumination. The ParFocus zoom system delivers a homogeneous, sharply-defined and reproducible laser spot on the retina, which minimizes heat-related side effects on the patient's cornea. The active ClearView physician safety filter offers not only a unique and true-to-color slit lamp image, but also reliably protects the physician, automatically swinging into position when the therapy beam is activated.

Navigate efficiently and confidently

The VISULAS 532s has a language-independent color touch screen for convenient and easy operation. The removable control panel can easily adjust to the viewing angle of the user. The PowerPress control enables direct selection of the power setting, without losing sight of the patient's fundus.



Navigating with fingertip precision:
the electronic micromanipulator.



... Removable control panel ...




... can be placed ergonomically
on the instrument table.

VISULAS 532s VITE

Expanding with your needs.
Fast. Flexible. Physician-controlled.

In addition to the conventional single-spot mode, the VISULAS® 532s VITE can also operate in multi-spot mode. A linear sequence of up to 12 laser pulses can be triggered at the touch of a button on the slit lamp joystick of the VISULAS 532s VITE.



3D illustration of a contoured VITE cascade



Fresh and pigmented laser lesions
(pulse duration: 20 ms)²⁾



Significantly reduce treatment times

The automated micromanipulator reliably controls the fast progression of an entire sequence of laser spots. Extremely precise, and long-term stable linear motors move the laser beam in just a few milliseconds to the next pre-calculated target position reducing conventional treatment time by 30 % to 60 %¹⁾.

Customize treatment strategies


The active control on the trigger button of the slit lamp allows the physician to maintain control over treatment progress at all times, with "fingertip precision". A slight rotation or movement of the joystick is all that is needed to adjust the position of the aiming beam as treatment proceeds. Flexible and precise. A multi-spot cascade allows the treatment strategy to be customized to the irregular contours of the retina.

Intelligently streamline treatment workflows

The clear layout of the graphic elements on the control panel follows a logical sequence and efficiently supports clinical workflows. Multi-spot treatment parameters, such as spot spacing, number of laser spots per sequence, and angle of orientation of a sequence, are homogeneously embedded in the VISULAS 532s user interface and can be adjusted during treatment. The illuminated SPOTview display allows the physician to continuously monitor the laser spot diameter, even in darkened environments.

Handle treatment interruptions with no stress

The physician may interrupt treatment at any time in an instant: with spot precision within the current sequence and in a controlled manner. This ensures the physician maximum precision and the patient optimum protection.



Comparison of average treatment time per 100 laser lesions for PRP treatments in conventional single-pulse mode of the VISULAS 532s (Group A: 35 patients) compared with PRP treatments in multi-spot mode of the VISULAS 532s VITE (Group B2: 37 patients). Time saving of around 40 %¹⁾.

¹⁾ Roeckl A, Blum M: Multispot laser photoocoagulation with the VISULAS 532s VITE:

A comparative study of 101 patients, Publication no. 000000-1839-880, LAS.2979, 06/2010.


²⁾ Roeckl A, Blum M: Pan-retinal laser photoocoagulation with reduced pulse duration – first experiences with linear spot cascades, Klinische Monatsblätter der Augenheilkunde, 2011 (DOI: 10.1055/s-0031-1273432)

VISULAS 532s VITE


Results to be proud of.

Effective. Gentle. Comfortable.

The VISULAS® 532s VITE offers a clinically effective and gentle retinal laser therapy for conventional treatment strategies, such as pan-retinal photo-coagulation, central grid coagulation and central focal coagulation.



3D illustration of a circular VITE cascade



Combined grid/PRP treatment
with VITE cascades³⁾

The must-have feature: conventional therapy with single pulses

The VISULAS 532s provides effective single-spot treatment using laser parameters that have established themselves as the evidence-based gold standard in accordance with the results of major clinical studies for the treatment of diabetic retinopathy and diabetic maculopathy (DRS, ETDRS, mETDRS, DRCR.net).


The exclusive add-on: conventional therapy with spot sequences

The particular advantage of the VISULAS 532s VITE is that it is the only coagulation laser that conforms to study-recommended laser settings, not only in the single-spot mode but also in the multi-spot mode (e.g. 50 ms pulse duration or 50 µm laser beam diameter). Initial clinical results prove that a significant reduction in treatment time can be achieved – without compromising the gold standard.


The compelling highlight: gentle therapy with short pulses

With the VISULAS 532s, it is also possible in short-pulse mode to set laser lesions efficiently and in a controlled fashion, with typical laser pulse duration between 10 and 40 ms – thanks to guaranteed laser output of 1.5 watt. Although a higher laser output is necessary to compensate for the lower pulse duration, in order to achieve a lesion of the same intensity, short-pulse photocoagulation is a gentle alternative to conventional laser therapy. This is due to a significantly lower energy load per unit per area².

In addition, shorter laser pulses cause less damage to surrounding tissue, since thermal conduction is not able to progress as far during the shorter pulse durations, thus protecting the inner layers of the retina (see images). This significantly enhances patient comfort²⁾. Short pulses create sharply outlined lesions with highly controlled heat propagation. Clinical research is currently underway on the potential positive effects of the combination of short pulse durations and small beam diameters.



OCT B scans (Cirrus™ HD-OCT by Carl Zeiss Meditec, Inc.) prove: A 20 ms laser pulse shows a less extensive lesion than a 200 ms pulse³⁾ (the beam diameter was 100 µm in both cases).



²⁾ Roeckl A, Blum M: Pan-retinal laser photoocoagulation with reduced pulse duration – first experiences with linear spot cascades, Klinische Monatsschriften der Augenheilkunde, 2011 (DOI: 10.1055/s-0031-1273432)

³⁾ Dr. Fang Lu, West China Hospital, Sichuan University, Chengdu, China

VISULAS 532s

A credible team player in operative use.
Adaptable. Portable. Versatile.

The VISULAS® 532s is designed for universal use. Due to its adaptable concept, the VISULAS 532s is ready to use in a blink of an eye – for the outpatient retina clinic or operating room. A wide range of applicators and high-quality accessories compliment the high-performance, multifaceted VISULAS 532s.

VISULINK 532/U – compatible with many diagnostic slit lamps

The VISULINK® 532/U – consisting of optical laser link and mechanical adapter – offers spot sizes of 50 µm to 1000 µm, and is equipped with a true-to-color physician safety filter. It converts diagnostic slit lamps made by Carl Zeiss, or e.g. Haag-Streit into fully-fledged laser workstations in an instant. The easy to use mechanism allows fast switching between different diagnostic slit lamps. To use the slit lamp diagnostically, the VISULINK 532/U can be swiveled out of the way.



VISULINK 532/U in operation

LIO 532s/Trion – high contrast and reliable

The Heine indirect ophthalmoscope specifically is modified for the application requirements of Carl Zeiss: the LIO 532s/Trion is suitable for both the VISULAS 532s and the multi-wavelength laser VISULAS Trion. It is light and stable and stands out in particular due to its high aiming beam contrast compared to the retinal background.



LIO 532s/Trion indirect ophthalmoscope

VISULAS 532s – strong performance in the operating room

With the coagulation laser VISULAS 532s and the surgical microscope OPMI LUMERA® 700 Carl Zeiss provides surgeons with a perfect duo which interact in harmony, guaranteeing successful surgical interventions in the posterior eye segment. The variable view, font size and color contrast of VISULAS 532s display are predestined for application in the operating room: all parameters are always clear and distinctly legible – even from a distance and in a darkened environment.



VISULAS 532s and OPMI LUMERA 700
in the operating room

The precise centration of the aiming and therapeutic beam, combined with the integrated slit illumination of the OPMI LUMERA 700, create an optimally illuminated treatment field and thus guarantees confident operation, whether in the periphery or in the vicinity of the macula. The fixed ClearView physician safety filter only minimally increases the height of the view for the operator. The active physician safety filter always provides effective protection when the laser is triggered. It gives the operator unfiltered vision when the laser is not in use.



Laser endoprobe with protective cap

VISULAS 532s

A credible team player in operative use.
Adaptable. Portable. Versatile.

The VISULAS® 532s is designed for universal use. Due to its adaptable concept, the VISULAS 532s is ready to use in a blink of an eye – for the outpatient retina clinic or operating room. A wide range of applicators and high-quality accessories compliment the high-performance, multifaceted VISULAS 532s.

VISULINK 532/U – compatible with many diagnostic slit lamps


The VISULINK® 532/U – consisting of optical laser link and mechanical adapter – offers spot sizes of 50 µm to 1000 µm, and is equipped with a true-to-color physician safety filter. It converts diagnostic slit lamps made by Carl Zeiss, or e.g. Haag-Streit into fully-fledged laser workstations in an instant. The easy to use mechanism allows fast switching between different diagnostic slit lamps. To use the slit lamp diagnostically, the VISULINK 532/U can be swiveled out of the way.

LIO 532s/Trion – high contrast and reliable

The Heine indirect ophthalmoscope specifically is modified for the application requirements of Carl Zeiss: the LIO 532s/Trion is suitable for both the VISULAS 532s and the multi-wavelength laser VISULAS Trion. It is light and stable and stands out in particular due to its high aiming beam contrast compared to the retinal background.



VISULINK 532/U in operation



LIO 532s/Trion indirect ophthalmoscope

Technical data

VISULAS 532s

VISULAS 532s / VISULAS 532s VITE	
Laser type	Frequency doubled solid state laser
Wavelength	532 nm
Aiming beam	Diode, 620 to 650 nm, adjustable brightness max. 1 mW at the cornea
Rated voltage and frequency	100 V to 240 V, 50/60 Hz
Pulse duration (single pulse)	10 to 2500 ms, cw
Pulse duration (VITE option)	20 to 50 ms
Pulse interval (single pulse)	100 to 6000 ms
Max. power	1.5 watts at the cornea
Cooling system	Thermoelectric
Laser console dimensions	H 150 mm x W 300 mm x D 400 mm (H 59 x W 118 x D 157) inches
Weight	14 kg (30.8 lbs)
Accessories	LIO 532s/Trion indirect ophthalmoscope, Instrument table, laser safety goggles, contact lenses, laser warning light

LSL 532s / LSL 532s VITE laser slit lamp	
Laser beam delivery	Interlaced with slit illumination system
Laser treatment spot size	continuously adjustable, 50–1,000 µm (without contact lens), parfocal, larger spot sizes depending on contact lens used
Illumination	12 V, 30 W brightness continuously adjustable
Slit adjustment	Slit length: variable in steps of 1/3/5/9/14 mm Slit width: continuously adjustable from 0 to 14 mm Slit image rotation: 0°, ±45°, 90°
Magnification	5 x, 8 x, 12 x, 20 x, 32 x
Physician safety filter	True-to-color, ClearView
Micromanipulator	Servo-electric
Weight	11 kg (24.2 lbs)
Accessories	Tonometer, co-observation tube, video documentation equipment from the range of accessories for the slit lamps SL120 and SL130

VISULINK 532/U	
Laser treatment spot size	continuously adjustable, 50–1,000 µm (without contact lens), parfocal
Compatible slit lamps	ZEISS SL 115 Classic, SL 120, SL 130, 20 SL, 30 SL Haag Streit 900® BM / 900® BQ
Physician safety filter	True-to-color
Weight	0.4 kg (0.9 lbs)
Accessories	Transport case for VISULINK 532/U and VISULAS 532s



Your local contact:

Argentina
 Carl Zeiss Argentina S.A.
 Calle Nahuel Huapi 4015 / 25
 C1430 BCO Buenos Aires
 Argentina
 Phone: +54 11 45 45 66 61
 bruzzi@zeiss.com.ar

Australia
 Carl Zeiss Pty. Ltd.
 Unit 13, 2 Eden Park Drive
 North Ryde, New South Wales 2113
 Australia
 Phone: +61 2 9020 1333
 med@zeiss.com.au

Austria
 Carl Zeiss GmbH
 Laxenburger Str. 2
 1100 Vienna
 Austria
 Phone: +43 1 79 51 80
 austria@zeiss.org

Belgium
 Carl Zeiss NV-SA
 Ikaroslaan 49
 1930 Zaventem
 Belgium
 Phone: +32 2 719 39 11
 info@zeiss.be

Brazil
 Carl Zeiss do Brasil Ltda.
 Av. Nações Unidas, 21711
 CEP04795-100 São Paulo
 Brazil
 Phone: +55 11 5693 5521
 medbrasil@zeiss.org

Canada
 Carl Zeiss Canada Ltd.
 45 Valleybrook Drive
 Toronto, ON M3B 2S6
 Canada
 Phone: +1 800 387 8037
 micro@zeiss.com

China
 Carl Zeiss Shanghai Co. Ltd.
 1/f., Ke Yuan Building
 11 Ri Yin Nan Road
 Waigaoqiao Free Trade Zone
 2005 Yang Gao Bei Road
 Shanghai 200131
 China
 Phone: +86 21 5048 17 17
 sro@zeiss.com.cn

Czech Republic
 Carl Zeiss spol. s.r.o.
 Radlická 14/3201
 150 00 Prague 5
 Czech Republic
 Phone: +420 233 101 221
 zeiss@zeiss.cz

France
 Carl Zeiss Meditec France SAS
 60, route de Sartrouville
 78230 Le Pecq
 France
 Phone: +33 1 34 80 21 00
 med@zeiss.fr

Germany
 Carl Zeiss Meditec VG mbH
 Carl-Zeiss-Strasse 22
 73446 Oberkochen
 Germany
 Phone: +49 7364 20 6000
 vertrieb@meditec.zeiss.com
 Surgical Ophthalmology:
 Phone: +49 800 470 50 30
 iol.order@meditec.zeiss.com

Hong Kong
 Carl Zeiss Far East Co. Ltd.
 Units 11-12. 25/F
 Tower 2, Ever Gain Plaza
 No. 88 Container Port Road
 Kwai Chung
 Hong Kong
 Phone: +852 2332 0402
 czfe@zeiss.com.hk

India
 Carl Zeiss India Pvt. Ltd.
 22. Kensington Road
 Ulsoor
 Bangalore 560 008
 India
 Phone: +91 80 2557 88 88
 info@zeiss.co.in

Italy
 Carl Zeiss S.p.A.
 Viale delle Industrie 20
 20020 Arese (Milan)
 Italy
 Phone: +39 02 93773 1
 post@zeiss.it

Japan
 Carl Zeiss Meditec Japan Co. Ltd.
 Shinjuku Ku
 Tokyo 160-0003
 22 Honchio-Cho
 Japan
 Ophthalmic instruments:
 Phone: +81 3 33 55 0331
 medsales@zeiss.co.jp
 Surgical instruments:
 Phone: +81 3 33 55 0341
 cmskoho@zeiss.co.jp

Malaysia
 Carl Zeiss Sdn Bhd.
 Lot2, Jalan 243/51 A
 46100 Petaling Jaya
 Selangor Darul Ehsan
 Malaysia
 Phone: +60 3 7877 50 58
 malaysia@zeiss.com.sg

Mexico
 Carl Zeiss de México S.A. de C.V.
 Avenida Miguel Angel de Quevedo
 496
 04010 Mexico City
 Mexico
 Phone: +52 55 59 99 0200
 cz-mexico@zeiss.org

Netherlands
 Carl Zeiss B.V.
 Trapezium 300
 Postbus 310
 3364 DL Sliedrecht
 Netherlands
 Phone: +31 184 43 34 00
 info@zeiss.nl

New Zealand
 Carl Zeiss (N.Z.) Ltd.
 15B Paramount Drive
 P.O. Box 121 - 1001
 Henderson, Auckland 0650
 New Zealand
 Phone: +64 9 838 5626
 med@zeiss.com.au

Poland
 Carl Zeiss sp. z o.o.
 ul. Lopuszanska 32
 02-220 Warsaw
 Poland
 Phone: +48 22 858 2343
 medycyna@zeiss.pl

Singapore
 Carl Zeiss Ltd.
 50 Kaki Bukit Place
 Singapore 415926
 Singapore
 Phone: +65 6741 9600
 info@zeiss.com.sg

South Africa
 Carl Zeiss (Pty.) Ltd.
 363 Oak Avenue
 Ferndale
 Randburg 2194
 South Africa
 Phone: +27 11 886 9510
 info@zeiss.co.za

South Korea
 Carl Zeiss Co. Ltd.
 Seoul 121-828
 Mapo-gu
 141-1, Sangsu-dong
 2F, BR Elitel Bldg.
 South Korea
 Phone: +82 2 3140 2600
 korea@zeiss.co.kr

Spain
 Carl Zeiss Meditec Iberia S.A.
 Ronda de Poniente, 15
 Tres Cantos
 28760 Madrid
 Spain
 Phone: +34 91 203 37 00
 info@zeiss.es

Sweden
 Carl Zeiss AB
 Tegeluddsvägen 76
 10254 Stockholm
 Sweden
 Phone: +46 84 59 25 00
 info@zeiss.se

Switzerland
 Carl Zeiss AG
 Feldbachstrasse 81
 8714 Feldbach
 Switzerland
 Phone: +41 55 254 7534
 med@zeiss.ch

Thailand
 Carl Zeiss Thailand
 Floor 8, Thosapol Land Building 2
 230 Ratchadapisek Road
 Huaykwang, Bangkok 10310
 Thailand
 Phone: +66 2 274 06 43
 thailand@zeiss.com.sg

United Kingdom
 Carl Zeiss Ltd.
 15-20 Woodfield Road
 Welwyn Garden City
 Hertfordshire, AL7 1JQ
 United Kingdom
 Phone: +44 1707 871200
 info@zeiss.co.uk

United States of America
 Carl Zeiss Meditec, Inc.
 5160 Hacienda Drive
 Dublin, CA 94568
 USA
 Phone: +1 925 557 4100
 info@meditec.zeiss.com



Manufacturer:

Carl Zeiss Meditec AG
 Goeschwitzer Strasse 51-52
 07745 Jena
 Germany
www.meditec.zeiss.com/visulas532s